

OOSS THz First Generation Light Wave Instrument Product Manual

1. Introduction:

- This product manual is intended to provide you with detailed information about the "THz Light Wave Instrument" (hereinafter referred to as the "Product") to ensure that you can use the product correctly and safely. Please read this manual carefully before use and follow all operating and safety guidelines.
- Please note that this product is an electrical appliance, improper use may cause damage or safety risks. Therefore, it is essential to understand and follow all instructions in this manual.

2. Product Overview:

- This product features advanced design and functionality, combining terahertz waves and thermal energy to provide convenient and efficient cellular therapy. Key functions include promoting cellular metabolism, enhancing cellular immunity, delaying cellular aging, and more.
- This product is not suitable for individuals with metal implants or using electromagnetic products internally or externally, those experiencing acute illnesses, open fractures caused by external force injuries, congenital heart valve defects, etc.
- Refer to the accompanying charts to understand the names and functions of the components of this product.

3. Technical Specifications:

- Voltage: 110-130V/220-250V
- Frequency: 50Hz/60Hz
- Power: 750W
- Dimensions: (Length) 347mm
- Gear: 1-3

4. Safety Instructions:

- Handle this product with care due to the presence of fragile crystals to avoid impact or dropping.
- Do not block the air inlet and outlet during use to prevent overheating.
- Avoid immersing the product in water or contact with other liquids to prevent electric shock.
- Remove all metal accessories from your body before using the product. Remove contact lenses before facial use.
- Ensure that the plug and socket are dry and inserted correctly during use.
- Avoid placing the product too close to any part of the body during use to prevent burns or other injuries.
- Unplug the power plug when the product is not in use.
- Do not use the product outdoors. The product is for household use only. Do not use the product for other purposes.
- Do not use the product for an extended period (**each therapy session should not exceed 15 minutes**). The product is equipped with a resettable thermal link to protect the motor from overheating. If the motor overheats, the link will disconnect the power and automatically reconnect it after the motor cools down.

- If the power cord is damaged, it must be replaced by the manufacturer, its service agent, or similar qualified persons to avoid danger.
- This product is not suitable for use by individuals (including children) with reduced physical, sensory, or mental capabilities or lack of experience and knowledge, unless they have been given supervision or instruction concerning the use of the device by a person responsible for their safety. Children should be supervised to ensure their safety.
- Do not leave the product unattended after turning it on.
- If the product cannot be turned on, it is in a malfunction state. Please contact the local service center.

5. Assembly and Setup:

- After opening the box, take out the two components.
- Before assembling the product, ensure that the power is disconnected, and check that all parts are intact.
- As shown in the diagram, align (1) the electronic tube with the part indicated by the arrow and install it onto (2) the base. Slightly rotate to adjust to ensure a firm connection. After both parts are locked, they will make a "click" sound.
- Insert the power plug and toggle the silver switch to adjust the gear and power size.



6. Operating Instructions and Usage Tips:

- Gear 1: High wind speed, low thermal energy; Gear 2: Low wind speed, medium thermal energy; Gear 3: High wind speed, high thermal energy.
- Drink a glass of water/terahertz water at the beginning and end of each session.
- Each use should last for 15 minutes. Multiple uses per day (3 to 5 times) are encouraged.
- During multiple uses, ensure that the product has completely cooled down after every 15 minutes before continuing use.
- During daytime use, blow from head to toe, blowing the internal organs to quickly replenish energy.
- During nighttime use, blow from the knees to the soles of the feet to calm the nerves and aid sleep.
- In painful areas/main organs/acupoints, blow each area for at least 5 minutes.
- After finishing, we also encourage you to conduct an experiment together. Please silently recite: "Please remove all negative energy from the body and pray to the OOSS Prayer Hall."

7. Maintenance and Cleaning:

- Regularly check the power cord and plug of this product to ensure that there is no wear or damage.
- Do not use chemical cleaners or solvents to clean this product and its parts to avoid surface damage.

8. Troubleshooting:

- If the electronic tube is broken or damaged, please contact the local service center for replacement.
- If (1) the electronic tube and (2) the base cannot be connected smoothly, check if the transparent crystal tube inside the electronic tube is misaligned.
- If there is a fault or sparks at the air inlet and wire connection during use, please contact the local service center for replacement.

9. Safety Certification and Compliance: (To be confirmed)

- Information about any safety certifications or standards the product meets.
- Compliance with relevant regulations and guidelines.

10. Warranty Information:

- The warranty period of this product is 1 year (from the date of purchase).
- During the warranty period, any failures caused by the product itself, please send it to the local service center for warranty service.
- After the warranty period expires, you will be responsible for the repair costs, and the company will provide repair and maintenance services.
- Warranty is void under the following circumstances:
 1. Exceeding the warranty period.
 2. Failure caused by not operating the device according to the instructions in the manual.
 3. The following situations are not within the scope of warranty:
 - a. Damage caused by improper use, storage, or maintenance by the consumer.
 - b. Damage caused by disassembly and repair by non-company designated repair departments.
 - c. Damage caused by fall and impact.

11. Other Information:

- Environmental protection period for this product:

The environmental protection period for this electronic appliance is 10 years. The environmental protection period refers to the period during which electronic appliances do not cause serious pollution to the environment or serious harm to human health or property.